
activforce2



Meet the handheld digital Dynamometer & Inclinometer in one

Powerful Technology in a Portable Package

The Activforce 2 is an advanced, patented digital strength and angle measurement platform that offers an innovative, portable, data-driven solution at an affordable price while increasing retention rates and driving better client outcomes





Affordable and Effective

Combining the functionality of a digital dynamometer and inclinometer in one, Activforce 2 enables physiotherapists to drive better outcomes by creating and prescribing recovery programs supported by accurate, conclusive data.



- Activforce 2 at an affordable price making dynamometry available to most physiotherapists.
- Generate accurate and conclusive data to increase client engagement and improve plan of care
- Measure force and range of motion with a built in inclinometer




Accurately Measure

- Active and Passive Range of Motion (AROM & PROM)
- Peak Force
- Average Force
- Force over Time
- Force/ Weight Ratio
- Muscle Symmetry



Quantify Progress

- Results in real time
- ROM progress
- Strength progress
- Data from test session reports delivered to your inbox



•

Average values

Test 1	19.38 kg
Test 2	19.64 kg

Peak Force (kg)

Right	23.78 kg
Left	21.95 kg
Strength Difference	1.84 kg
Percentage Difference	8.03%

Shoulder extension

AROM Symmetry	
Right	224.08°
Left	205.08°
Angle Difference	19.01°
Percentage Difference	8.86%

Proven Results

One of the biggest impediments to a successful treatment is client abandonment before completion of a prescribed course of care. Activforce 2 can help improve compliance by offering a more accurate and objective gauge of the effectiveness of care.

arli ee Wil	FiCall 🗢 12: Sumi	04	
S	Symmetry	All Data	
Right Hand grip			
1	20.19 kg Peak force	17.32 kg Avg force	
2	19.99 kg Peak force	13.80 kg Avg force	
Left H	and grip		
1	18.58 kg Peak force	16.73 kg Avg force	
2	17.79 kg Peak force	16.14 kg Avg force	
	Email Su	immary	
F	Perform Anothe	r Measurement	





Increase patient engagement and retention with objective, quantifiable data.

The new **Activforce 2 device + companion app** allows you to perform both strength and range of motion tests so you can measure and track your patients' progress.

- Small, sleek and portable (4.6 ounces, 3.7 inches long to be exact)
- Durable
- Measures up to 200 lbs in strength
- Force accuracy within +/- 1.4 lb + 5% of the applied force
- Bluetooth-enabled, no internet required
- No charging necessary, comes with AA batteries that will last 6-12 months depending on usage

What's Included:



- New Activforce 2 device with both strength and angle measurement capabilities
- High quality carrying case that comes complete with:
 - New locking mechanism which allows you to easily switch attachments that better contour to the patient's body part. This locking mechanism will also allow the PT to purchase additional attachments in the future without having to buy a new device.
 - Rounded pad for testing most body parts
 - Arched pad to comfortably test angular body parts
 - Hand and leg straps to hold the device securely
 - A belt to aid in performing pull measurements
 - A belt loop plate to hold the straps and belt securely to the device
- A free, lifetime basic account
- Complementary knowledge base including strength testing guides and how-to videos

Accessories





One Comprehensive Portable Kit

Strength Testing Protocol Samples



Flexion

Patient:Seated or SupineShoulder:Flexed to 90 degreesDevice:Proximal to elbow



Abduction

 Patient:
 Prone

 Shoulder:
 Abducted to 90 degrees

 Elbow:
 Flexed to 90 degrees

 Device:
 Base of head



Internal Rotation

 Patient:
 Prone

 Shoulder:
 Abducted to 90 degrees, externally rotate to comfortable position

 Elbow:
 Flexed to 90 degrees

 Device:
 Proximal to wrist



Knee Extension

Patient:SeatedKnee:Extended to 0 degreesDevice:Mid-shin



Hip Flexion

Patient:	Prone
Hip:	Flexed to 100 degrees
Knee:	Flexed to 90 degrees
Device:	Mid-thigh

Flexion

 Patient Position:
 Seated

 Neck Position:
 Flexed to comfortable angle

 Device Position:
 Forehead



Rotation

Patient Position:SeatedNeck Position:Slight rotationDevice Position:Temple

Appendix

Introduction to Activforce 2 Basic Assessments

Activforce 2 Strength Measurement

Activforce 2 Attachment

Activforce 2 Hand Strap

Activforce 2 Belt Demonstration

Activforce 2 Range of Motion Measurement