



HYPERVOLT

Featuring Quiet Glide™ Technology

NEW FEATURE

Now with pressure sensor technology

The world's most advanced percussion massage device just got better. The Hypervolt now features three LED lights to indicate the pressure level for immediate feedback for your practitioner / massage therapist.



ACTIVATING THE PRESSURE SENSOR

1. Power ON (bottom of handle)
2. Press "speed setting button" (level 1) and wait 5 seconds to activate the pressure sensor
3. Select desired speed setting (press once per speed)

Pressure Sensor LED Lights (3)

HEAD ATTACHMENTS



Fork



Bullet



Flat



Cushion



Round