


Active compression squeezes greater benefits out of cold therapy.

Infrared thermal imagery shows that Game Ready's patented active compression technology enhances the delivery of cold therapy.

Head to head and knee to knee, Game Ready® provides faster, deeper, and longer-lasting cold than other cryotherapy products.

Game Ready takes injury recovery to a higher level by combining active pneumatic compression and circumferential cold therapy in one convenient system. Our active compression not only provides greater therapeutic benefits over cold alone, but shapes our specialized wraps to the contours of the body—increasing surface contact, helping to optimize cold therapy delivery.

More effective cold. More effective compression. Game Ready.



Side view shown. First skin temperature measurement immediately after completing 30 minutes of treatment, subsequent images at 10 and 30 minutes post-treatment. All devices applied per manufacturer's instructions and all measurements taken from uniform distance. All devices applied on different days allowing full recovery for test subject.

Active Compression and Cold Therapy for Greater Therapeutic Benefits

Game Ready simultaneously delivers active (intermittent) compression and adjustable cold therapy for faster, deeper, longer-lasting cooling (see *thermal imagery*).

Active compression helps facilitate the healing process, reduces swelling, **and** conforms ergonomic, flexible wraps to body contours for increased surface contact and improved heat/cold transfer.

Patented circumferential wraps fully surround injury site, covering more surface area, further enhancing the delivery of cold therapy.

Game Ready Wraps reach key patient pain points for ACL repair or total knee replacement such as the upper shin and back of the knee, respectively.